

DG Gel[®] 8 Cards

Unique 8-column design offering clear, easy-to-read results

Column Agglutination Method Using a Gel Matrix



Proven performance

- High sensitivity and specificity
- Gel agglutination technique provides proven accuracy and reliability
- Consistent and reproducible interpretation of results
- Reactions remain stable and allow for additional review/QA

Customizable workflow

- Multiple test configurations
- For use in automated, semi-automated and manual testing

Test configurations


- ABO/Rh (D) blood group typing (forward and reverse)
- Antibody screening
- Antibody identification
- Crossmatch: immediate spin and IgG
- Direct antiglobulin technique
- Rh phenotyping
- Kell typing
- Donor confirmation
- Weak D

TYPING

DG Gel® 8 Cards

Column Agglutination Method Using a Gel Matrix

Technical specifications



DG Gel® 8 A/B/D (x2)

Configuration: A/B/D/Ctl/A/B/D/Ctl
Presentation: 2x25 cards

Reference # 210384



DG Gel® 8 Neutral

Configuration: N/N/N/N/N/N/N/N
Presentation: 2x25 cards


Reference # 210393



DG Gel® 8 Anti-IgG (Rabbit)

Configuration: IgG/IgG/IgG/IgG/IgG/IgG/IgG/IgG
Presentation: 2x25 cards


Reference # 210394



DG Gel® 8 ABO/Rh 2D

Configuration: A/B/A,B/DVI-/DVI+/Ctl/N/N
Presentation: 2x25 cards

Reference # 210128



DG Gel® 8 AB (x4)

Configuration: A/B/A/B/A/B/A/B
Presentation: 2x25 cards


Reference # 210386



DG Gel® 8 Anti-AB (x8)

Configuration: AB/AB/AB/AB/AB/AB/AB/AB
Presentation: 1x25 cards

Reference # 210200



DG Gel® 8 Direct Coombs (x2)

Configuration: Poly/IgG/C3d/Ctrl/Poly/IgG/C3d/Ctrl
Presentation: 1x25 cards

Reference # 210129



DG Gel® 8 Anti-K (x8)

Configuration: K/K/K/K/K/K/K/K
Presentation: 1x25 cards


Reference # 210225



DG Gel® 8 Anti-D (x8)

Configuration: D/D/D/D/D/D/D/D
Presentation: 1x25 cards


Reference # 210202



DG Gel® 8 Rh Pheno

Configuration: C/E/c/e/C/E/c/e
Presentation: 1x25 cards

Reference # 210228



DG Gel® 8 Anti-C (x8)

Configuration: C/C/C/C/C/C/C/C
Presentation: 1x25 cards


Reference # 210201



DG Gel® 8 Anti-E (x8)

Configuration: E/E/E/E/E/E/E/E
Presentation: 1x25 cards


Reference # 210224



DG Gel® 8 Anti-c (x8)

Configuration: c/c/c/c/c/c/c/c
Presentation: 1x25 cards

Reference # 210226



DG Gel® 8 Anti-e (x8)

Configuration: e/e/e/e/e/e/e/e
Presentation: 1x25 cards

Reference # 210227